Black

Date	27/04/2020
Client	MOLNLYCKE
Description	IFU MELGISORB Ag EU
Dressing size	5X5, 10X10, 15X15, 20X30 3X44
Client Ref	256055, 256105, 256155, 256205, 256605
Finished good	10015689 to 5693
No of Colours	1
Supplier details	
Item number	10015704 - 21/18/E4
Size flat	95mm (W) x 139mm (h)
Final Folded Size	60 PAGE BOOKLET
View when folded	MELGISORB AG
Print process	Offset Litho
Raw Material ref	RMS006 -80gsm
Version	A DBT


 AMS Studio, Advanced Medical Solutions Ltd.,
Road One, 33 Premier Park, Winsford Industrial Estate,
Winsford, Cheshire CW7 3RT E: ams.studio@admedsol.com
T: +44 (0) 1606 54 5520 F: +44 (0) 1606 54 5626

SUPPLIER

Please acknowledge receipt of artwork as this legend will be the final specification. Check details and sign / stamp if this legend detail and artwork is correct and can be printed on the substrate indicated. Proceed to repro and proof.

Do not print until you receive an official order
Do not destroy any previous plates until
AMS confirm it is ok for you to do so. Any
plates destroyed without this confirmation
will be replaced at the supplier's cost.

Melgisorb® Ag

STERILE R

 Single use
Usage unique

CE 2797

w/w 0.5% Ag

Distributed by:
Mölnlycke Health Care AB

Manufacturer
 **Advanced Medical Solutions Ltd**
Premier Park, 33 Road One,
Winsford Industrial Estate, Winsford, Cheshire, CW7 3RT, UK

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

Product Description

Melgisorb Ag is a highly absorbent, sterile, non woven, antimicrobial pad composed of a high G (guluronic acid) calcium alginate, carboxymethylcellulose (CMC), and an ionic silver complex (Silver Sodium Hydrogen Zirconium Phosphate, w/w 0.5% Ag) which releases silver ions in the presence of wound exudate. As wound exudate is absorbed, the alginate forms a gel, which aids autolytic debridement whilst maintaining a moist environment for optimal wound healing. The gel ensures intact removal without damage to healing tissue.

The silver ions released in the presence of wound exudate are an effective antimicrobial agent against a broad spectrum of microorganisms frequently associated with the bacterial colonisation and infection of wounds, for up to fourteen (14) days, based on in-vitro testing. Melgisorb Ag can also reduce the odour of infected wounds as a result of its effective antimicrobial properties.

Melgisorb Ag is an effective barrier to penetration by microorganisms.

Indications

Melgisorb Ag is indicated for use in the management of all moderate to heavily exuding partial to full thickness wounds including:

- post-operative wounds
- trauma wounds (dermal lesions, trauma injuries or incisions)
- leg ulcers
- pressure ulcers
- diabetic ulcers
- graft and donor sites
- cavity wounds
- superficial and partial thickness burns

As the Melgisorb Ag contains alginate, it may assist in supporting the control of minor bleeding in superficial wounds.

This product can also be used under compression bandages.

Melgisorb Ag may be used in the management of infected wounds or wounds in which there is an increased risk of infection under the medical supervision of a healthcare professional.

Melgisorb Ag is indicated for external use only

Contra-indications

Melgisorb Ag is not indicated for use on the following:

- dry or lightly exuding wounds
- individuals with a known sensitivity to alginates or silver
- to control heavy bleeding

Precautions

The dressing may adhere if used on dry or very lightly exuding wounds. If the dressing does adhere and is not easily removed, moisten with sterile saline solution prior to removal.

The dressing performance may be impaired by excess use of petroleum-based ointments.

Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).

The dressing must be removed prior to patients undergoing Magnetic Resonance Imaging (MRI) examinations.

In the event of clinical infection, topical silver does not replace the need for systemic therapy or other adequate infection treatment.

Clinicians / Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.

Melgisorb Ag is intended for short-term use, up to 30 days.

DIRECTIONS FOR USE

Site Preparation

- Debride when necessary and irrigate the wound site in accordance with standard protocols.
- Remove excess solution from surrounding skin.

Dressing Selection

- Select a size of Melgisorb Ag that is slightly larger than the wound.

Dressing Application

1. Cut (using sterile scissors) or fold the dressing to fit the wound. Discard any remaining dressing material due to the risk of contamination.
2. Apply to wound bed directly. Loosely fill deep wounds, ensuring the dressing does not overlap the wound margins.
3. Cover and secure Melgisorb Ag with a non-occlusive secondary dressing.

Dressing Change and Removal

1. Melgisorb Ag can remain in situ up to 7 days. Dressing change frequency will depend on wound condition and the level of exudate. Initially it may be necessary to change the dressing every 24 hours.
2. Reapply Melgisorb Ag when the secondary dressing has reached its absorbent capacity or whenever good wound care practice dictates that the dressing should be changed.
3. To reapply, gently remove the secondary dressing.
4. If the wound appears dry, saturate the dressing with normal sterile saline solution prior to removal.
5. Gently remove the dressing from the wound bed and dispose of according to local procedures and guidelines.
6. If required, irrigate the wound site in accordance with standard protocols prior to application of a new dressing.

Infected Wounds

1. If the wound is infected, appropriate medical treatment should be initiated.
2. Application of the Melgisorb Ag may be continued at the discretion of the attending medical practitioner/clinician.

Melgisorb Ag does not utilise animal tissue. Melgisorb Ag is not made with natural rubber latex.

Storage

Store at ambient temperature (below 25°C/ 77°F).

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CANADA Mölnlycke Health Care Inc.,
2010 Winston Park Drive, Suite 100,
Oakville, Ontario, L6H 5R7
Tel: +1 905 829 1502

CHINA Mölnlycke Health Care
(Shanghai) Co., Ltd.
RM 629, No.8 Hua Jing Road,
Wai Gao Qiao FTZ Shanghai, China
Tel: +86 10 5128 8571

CZECH REPUBLIC
Mölnlycke Health Care s.r.o.
Hajkova 2747/22, 130 00 Prague 3
Tel: +420 221 890 517 (Reception)
Tel: +420 221 890 511 (Customer Service)

DENMARK Mölnlycke Health Care ApS,
Gydevang 33, DK-3450 Allerød
Tel: +45 48 16 82 68
Customer Service Tel: +80 886 810

ESTONIA Mölnlycke Health Care OÜ,
Vabaõhumuuseumitee 4-67, EE-13522 Tallin
Tel: +372 671 1520

FINLAND Mölnlycke Health Care Oy,
Pitäjänmäentie 14, 3 krs, FI-00380 Helsinki
Tel: +358 201 622 300
Customer Service Tel: 0800 113 442

FRANCE Mölnlycke Health Care SAS,
13 Allée du Château Blanc
59290 WASQUEHAL
Tel: +33 320 122 555
Customer Service Tel: 0800 910 292

GERMANY Mölnlycke Health Care GmbH
Grafenberger Allee 297
DE-40237 Düsseldorf
Tel: +49 211 92 08 80
Customer Service Tel: 0800 186 21 80

HUNGARY Mölnlycke Health Care Kft
HU-1134 Budapest, Dévai u. 26-28. III. em.
Tel: +36 1 477 3080

INDIA Mölnlycke Health Care India Pvt Ltd
Universal Business Centre
Second Floor, Universal Trade Tower
Sector 49, Gurgaon Sohna Road
Gurgaon 122 018, Haryana
Tel: +91 124 4696 008

ITALY Mölnlycke Health Care s.r.l.,
Via Marsala 40/C, IT-21013 Gallarate (VA)
Tel: +39 0331 714 011

JAPAN Mölnlycke Health Care KK,
4F Concieria Nishi-Shinjuku Tower's West
6-20-7 Nishi-Shinjuku,
Shinjuku-Ku, Tokyo 160-0023
Tel: +81 3 6914 5004

KOREA Mölnlycke Health Care Korea Co., Ltd
2nd floor, Hansol Richville.
46 Wiryeseong-daero,
Songpa-gu, Seoul, 05627
South Korea
Tel: +82 2 3789 1402

LATVIA Mölnlycke Health Care AB,
Dzelzavas 120 M, Riga LV-1021
Tel: +372 6711 522

LITHUANIA Mölnlycke Health Care AB,
J. Galvydzio g. 3, LT 08236 Vilnius
Tel: +370 5 2274 58 37

NETHERLANDS Mölnlycke Health Care B.V.
Postbus 3196, 4800 DD Breda
Tel: +31 76 521 9663

NORWAY Mölnlycke Health Care AS
Postboks 6229 Etterstad, NO-0603 Oslo
Tel: +47 22 70 63 70
Customer Service Tel: 800 161 78

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SWITZERLAND Mölnlycke Health Care AG

Brandstrasse 24
CH-8952 Schlieren
Tel: +41 44 744 54 00
Customer Service
Tel German: 0800 563 195
Tel French: 0800 563 196

UNITED ARAB EMIRATES

Mölnlycke Health Care
Unit No. 603
Jumeirah Business Center 1, Plot No. G2
Jumeirah Lakes Towers, Dubai
Tel: +971 45 54 73 17

UNITED KINGDOM

Mölnlycke Health Care Ltd.
401 Grafton Gate
Milton Keynes
MK9 1AT
Tel: +44 870 60 60 766
Customer Service, Private
Tel: 0800 917 4919
Customer Service, Public sector
Tel: 0800 917 4918

USA Mölnlycke Health Care LLC

5550 Peachtree Parkway
Suite 500, Norcross, GA 30092
Tel: + 1-800-882-4582

Symbols used on labelling

Sterile
Method of sterilisation: Irradiation

Single use
Usage unique

Caution

Do not use if package is damaged
Ne pas utiliser si l'emballage est endommagé

Do not resterilize

Consult instructions
for use

Product code

Manufacturer

Batch number

Use by:
Year and Month

< 25°C
77°F

Store at ambient temperature

CE 2797

Authorized representative
in the European Community